

Anti-OGT Recombinant Rabbit Monoclonal Antibody



Catalog #: 1325

Aliases

OGT; O-Linked N-Acetylglucosamine (GlcNAc) Transferase; O-GLCNAC; HRNT1; OGT1; O-Linked N-Acetylglucosamine (GlcNAc) Transferase (UDP-N-Acetylglucosamine:Polypeptide-N-Acetylglucosaminyl Transferase); UDP-N-Acetylglucosamine--Peptide N-Acetylglucosaminyltransferase 110 KDa Subunit; UDP-N-Acetylglucosamine:Polypeptide-N-Acetylglucosaminyl Transferase; O-Linked N-Acetylglucosamine Transferase 110 KDa Subunit; O-GlcNAc Transferase Subunit P110; MGC22921; FLJ23071; Uridinediphospho-N-Acetylglucosamine:Polypeptide Beta-N-Acetylglucosaminyl Transferase; O-GlcNAc Transferase P110 Subunit; EC 2.4.1.255; EC 2.4.1.186; HINCUT-1; EC 2.4.1; XLID106; MRX106

Background

Gene Name: OGT

NCBI Gene Entry: [8473](#)

UniProt Entry: [O15294](#)

Application Information

Molecular Weight: Predicted, 117 kDa; observed, 110 kDa

Clonality: Rabbit monoclonal antibody

Clone ID: 23GB2505

Species Reactivity: Human, mouse, rat

Applications Tested: Western blotting (WB), flow cytometry (FCM), immunocytochemistry (IC)

Immunogen

A synthesized peptide derived from human OGT/O-Linked N-Acetylglucosamine Transferase

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:1,000-1:5,000

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

Anti-OGT Recombinant Rabbit Monoclonal Antibody


PAGE 2

Flow Cytometry (FCM): 1:2,000

Immunocytochemistry (IC): 1:100-1:1,000


Note: This product is for research use only.

Validation Data



Copyright ©2024 Genuin Biotechnologies LLC

Flow cytometric analysis of OGT expression in C2C12 cells using OGT antibody (Cat#1325, 1:2,000). Green, isotype control; red, OGT.



Copyright ©2024 Genuin Biotechnologies LLC

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS


SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

Anti-OGT Recombinant Rabbit Monoclonal Antibody

PAGE 3

Immunocytochemical staining of C2C12 cells with OGT antibody (Cat#1325, 1:1,000). Nuclei were stained blue with DAPI; OGT was stained magenta with Alexa Fluor® 647. Images were taken using Leica stellaris 5. Protein abundance based on laser Intensity and smart gain: Medium. Scale bar: 20 μ m.



Copyright ©2024 Genuin Biotechnologies LLC

Western blotting analysis using anti-OGT antibody (Cat#1325). Total cell lysates (30 μ g) from various cell lines were loaded and separated by SDS-PAGE. The blot was incubated with anti-OGT antibody (Cat#1325, 1:5,000) and HRP-conjugated goat anti-rabbit secondary antibody (Cat#201, 1:20,000) respectively. Image was developed using FeQ™ ECL Substrate Kit (Cat#226).

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM