

Anti-AMHR2 Recombinant Rabbit Monoclonal Antibody



Catalog #: 1420

Aliases

AMHR2; Anti-Mullerian Hormone Receptor Type 2; MISRII; MISR2; Muellerian Inhibiting Substance Type II Receptor; Anti-Muellerian Hormone Type II Receptor; Anti-Muellerian Hormone Type-2 Receptor; AMH Type II Receptor; MIS Type II Receptor; EC 2.7.11.30; AMHR; MRII; Mullerian; Inhibiting Substance Type II Receptor; Anti-Mullerian Hormone Receptor, Type II; Anti-Mullerian Hormone Receptor Type II

Background

Gene Name: AMHR2
NCBI Gene Entry: [269](#)
UniProt Entry: [Q16671](#)

Application Information

Molecular Weight: Predicted, 63 kDa; observed, 63 kDa
Clonality: Rabbit monoclonal antibody
Clone ID: 23GB3205
Species Reactivity: Rat
Applications Tested: Western blotting (WB), flow cytometry (FCM)

Immunogen

A synthesized peptide derived from human AMHR2

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:1,000-1:5,000
Flow Cytometry (FCM): 1:2,000

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041


WWW.GENUINBIOTECH.COM

Anti-AMHR2 Recombinant Rabbit Monoclonal Antibody

PAGE 2

Note: This product is for research use only.

Validation Data



Copyright ©2024 Genuin Biotechnologies LLC

Flow cytometric analysis of AMHR2 expression in H9c2 cells using AMHR2 antibody (Cat#1420, 1:2,000). Green, isotype control; red, AMHR2.


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Western blotting analysis using anti-AMHR2 antibody (Cat#1420). Total cell lysates (30 µg) from various cell lines were loaded and separated by SDS-PAGE. The blot was incubated with anti-AMHR2 antibody (Cat#1420, 1:5,000) and HRP-conjugated goat anti-rabbit secondary antibody (Cat#201, 1:20,000) respectively. Image was developed using FeQ™ ECL Substrate Kit (Cat#226).