

Anti-DDX17 Recombinant Rabbit Monoclonal Antibody



Catalog #: 1422

Aliases

DDX17; DEAD-Box Helicase 17; P72; DEAD/H (Asp-Glu-Ala-Asp/His) Box Polypeptide 17 (72kD); DEAD (Asp-Glu-Ala-Asp) Box Polypeptide 17; Probable ATP-Dependent RNA Helicase DDX17; DEAD (Asp-Glu-Ala-Asp) Box Helicase 17; RNA-Dependent Helicase P72; DEAD Box Protein P72; DEAD Box Protein P82; DEAD Box Protein 17; EC 3.6.4.13; EC 3.6.1; RH70

Background

Gene Name: DDX17

NCBI Gene Entry: [10521](#)

UniProt Entry: [Q92841](#)

Application Information

Molecular Weight: Predicted, 80 kDa; observed, 70 kDa

Clonality: Rabbit monoclonal antibody

Clone ID: 23GB3560

Species Reactivity: Human, mouse, rat

Applications Tested: Western blotting (WB), flow cytometry (FCM), immunocytochemistry (IC)

Immunogen

A synthesized peptide derived from human DDX17

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:1,000-1:5,000

Flow Cytometry (FCM): 1:2,000

Immunocytochemistry (IC): 1:100-1:1,000

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041



WWW.GENUINBIOTECH.COM

Anti-DDX17 Recombinant Rabbit Monoclonal Antibody

PAGE 2

Note: This product is for research use only.

Validation Data



SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS


SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

Anti-DDX17 Recombinant Rabbit Monoclonal Antibody

PAGE 3

Scale bar: 20 μ m.



Western blotting analysis using anti-DDX17 antibody (Cat#1422). Total cell lysates (30 μ g) from various cell lines were loaded and separated by SDS-PAGE. The blot was incubated with anti-DDX17 antibody (Cat#1422, 1:5,000) and HRP-conjugated goat anti-rabbit secondary antibody (Cat#201, 1:20,000) respectively. Image was developed using FeQ™ ECL Substrate Kit (Cat#226).

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM