

# Anti-Phospho-c-Jun (S73)+JunD (S100) Recombinant Rabbit Monoclonal Antibody



**Catalog #: 1794**

## Aliases

JUN; Jun Proto-Oncogene, AP-1 Transcription Factor Subunit; V-Jun Avian Sarcoma Virus 17 Oncogene Homolog; AP-1; Transcription Factor AP-1 Subunit Jun; Transcription Factor Jun; Proto-Oncogene C-Jun; Activator Protein 1; Jun Oncogene; C-Jun; AP1; P39; V-Jun Sarcoma Virus 17 Oncogene Homolog; Jun Activation Domain Binding Protein; Enhancer-Binding Protein AP1; Transcription Factor AP-1; Proto-Oncogene CJun; C-JUN; CJUN

## Background

Gene Name: JUN /JUND

NCBI Gene Entry: [3725](#)

UniProt Entry: [P05412/P17535](#)

## Application Information

Molecular Weight: Predicted, 36 kDa; observed, 38 kDa

Clonality: Rabbit monoclonal antibody

Clone ID: 23GB4495

Species Reactivity: Human, mouse, rat

Applications Tested: Western blotting (WB), flow cytometry (FCM), immunocytochemistry (IC)

## Immunogen

A synthesized peptide derived from human Phospho-c-Jun (S73)+JunD (S100)

## Isotype

Rabbit IgG

## Storage Buffer

Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.

## Storage

Store at -20 °C for one year.

## Recommended Dilutions

Western Blotting (WB): 1:1,000-1:5,000

Flow Cytometry (FCM): 1:2,000

Immunocytochemistry (IC): 1:100-1:1,000

---

### SUPPORT

SUPPORT@GENUINBIOTECH.COM  
TEL: +1-540-855-7041

### ORDERS

SALES@GENUINBIOTECH.COM  
FAX: +1-540-855-7041


[WWW.GENUINBIOTECH.COM](http://WWW.GENUINBIOTECH.COM)

## Anti-Phospho-c-Jun (S73)+JunD (S100) Recombinant Rabbit Monoclonal Antibody


PAGE 2

**Note:** This product is for research use only.

### Validation Data



Flow cytometric analysis of Phospho-c-Jun (S73)+JunD (S100) expression in C2C12 cells using Phospho-c-Jun (S73)+JunD (S100) antibody (Cat#1794, 1:2000). Green, isotype control; red, Phospho-c-Jun (S73)+JunD (S100).



Immunocytochemical staining of C2C12 cells with Phospho-c-Jun (S73)+JunD (S100) antibody (Cat#1794, 1:1,000). Nuclei were stained blue with DAPI; Phospho-c-Jun (S73)+JunD (S100) was

#### SUPPORT

SUPPORT@GENUINBIOTECH.COM  
TEL: +1-540-855-7041

#### ORDERS


SALES@GENUINBIOTECH.COM  
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

## Anti-Phospho-c-Jun (S73)+JunD (S100) Recombinant Rabbit Monoclonal Antibody

PAGE 3

stained magenta with Alexa Fluor® 647. Images were taken using Leica stellaris 5. Protein abundance based on laser Intensity and smart gain: Medium. Scale bar: 20 µm.



Western blotting analysis using anti-Phospho-c-Jun (S73)+JunD (S100) antibody (Cat#1794). Total cell lysates (30 µg) from various cell lines were loaded and separated by SDS-PAGE. The blot was incubated with anti-Phospho-c-Jun (S73)+JunD (S100) antibody (Cat#1794, 1:5,000) and HRP-conjugated goat anti-rabbit secondary antibody (Cat#201, 1:20,000) respectively. Image was developed using FeQ™ ECL Substrate Kit (Cat#226).

### SUPPORT

SUPPORT@GENUINBIOTECH.COM  
TEL: +1-540-855-7041

### ORDERS

SALES@GENUINBIOTECH.COM  
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM