

Anti-Histone deacetylase 10 Recombinant Rabbit Monoclonal Antibody



Catalog #: 2135

Aliases

HDAC10; Histone Deacetylase 10; Polyamine Deacetylase HDAC10; DKFZP761B039; HD10; EC 3.5.1.; EC 3.5.1.

Background

Gene Name: HDAC10

NCBI Gene Entry: [83933](#)

UniProt Entry: [Q969S8](#)

Application Information

Molecular Weight: Predicted, 71 kDa; observed, 75 kDa

Clonality: Rabbit monoclonal antibody

Clone ID: 23GB6165

Species Reactivity: Human, mouse, rat

Applications Tested: Western blotting (WB), flow cytometry (FCM), immunocytochemistry (IC)

Immunogen

A synthesized peptide derived from human HDAC10

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:1,000-1:5,000

Flow Cytometry (FCM): 1:2,000

Immunocytochemistry (IC): 1:100-1:1,000

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS


SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM




Copyright ©2024 Genuin Biotechnologies LLC

Flow cytometric analysis of Histone deacetylase 10 expression in HepG2 cells using Histone deacetylase 10 antibody (Cat#2135, 1:2,000). Green, isotype control; red, Histone deacetylase 10.



Immunocytochemical staining of HepG2 cells with anti-Histone deacetylase 10 antibody (Cat#2135, 1:1,000). Nuclei were stained blue with DAPI; Histone deacetylase 10 was stained magenta with Alexa Fluor® 647. Images were taken using Leica stellaris 5. Protein abundance based on laser Intensity and smart gain: Medium. Scale bar, 20 μ m.



Western blotting analysis using anti-Histone deacetylase 10 antibody (Cat#2135). Total cell lysates (30 µg) from various cell lines were loaded and separated by SDS-PAGE. The blot was incubated with anti-Histone deacetylase 10 antibody (Cat#2135, 1:5,000) and HRP-conjugated goat anti-rabbit secondary antibody (Cat#201, 1:20,000) respectively. Image was developed using FeQ™ ECL Substrate Kit (Cat#226).