

Anti-Homeobox A5 Recombinant Rabbit Monoclonal Antibody



Catalog #: 2180

Aliases

HOXA5; Homeobox A5; HOX1C; Homeobox Protein Hox-A5; Homeobox Protein Hox-1C; Homeo Box A5; HOX1; Homeobox Protein HOXA5; Homeo Box 1C; HOX1.3

Background

Gene Name: HOXA5

NCBI Gene Entry: [3202](#)

UniProt Entry: [P20719](#)

Application Information

Molecular Weight: Predicted, 29 kDa; observed, 43 kDa

Clonality: Rabbit monoclonal antibody

Clone ID: 23GB6380

Species Reactivity: Human

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human HOXA5

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:1,000-1:5,000

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Western blotting analysis using anti-Homeobox A5 antibody (Cat#2180). Total cell lysates (30 µg) from various cell lines were loaded and separated by SDS-PAGE. The blot was incubated with anti-Homeobox A5 antibody (Cat#2180, 1:5,000) and HRP-conjugated goat anti-rabbit secondary antibody (Cat#201, 1:20,000) respectively. Image was developed using FeQ™ ECL Substrate Kit (Cat#226).