

AF680-Conjugated Goat Anti-Mouse IgG (H&L) Secondary Antibody



Catalog #: 231

Application Information

Clonality: Goat Polyclonal antibody

Species Reactivity: Mouse

Applications Tested: ELISA (E), Western Blotting (WB), Immunohistochemistry (IHC),
Immunocytochemistry (IC), Flow Cytometry (FCM)

Conjugation: AF680

Immunogen

Mouse IgG

Storage Buffer

Supplied in PBS (pH 7.2) containing 50% glycerol, 1% BSA and 0.02% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:5,000-1:20,000

Immunohistochemistry (IHC): 1:100-1:1,000

Immunocytochemistry (IC): 1:100-1:1,000

Flow Cytometry (FCM): 1:1,000-1:4,000

Validation Data

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

**AF680-Conjugated Goat Anti-Mouse IgG
(H&L) Secondary Antibody**

PAGE 2



SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM