

Anti-Exportin 5 Recombinant Rabbit Monoclonal Antibody



Catalog #: 2741

Aliases

XPO5; Exportin 5; KIAA1291; Ran-Binding Protein 21; Exportin-5; Exp5; RANBP21

Background

Gene Name: XPO5

NCBI Gene Entry: [57510](#)

UniProt Entry: [Q9HAV4](#)

Application Information

Molecular Weight: Predicted, 136 kDa; observed, 136 kDa

Clonality: Rabbit monoclonal antibody

Clone ID: 24GB2610

Species Reactivity: Human, mouse, rat

Applications Tested: Western blotting (WB), flow cytometry (FCM), immunocytochemistry (IC)

Immunogen

A synthesized peptide derived from human Exportin 5

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:1,000-1:5,000

Flow Cytometry (FCM): 1:2,000

Immunocytochemistry (IC): 1:100-1:1,000

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041


ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041


WWW.GENUINBIOTECH.COM



Copyright ©2024 Genuin Biotechnologies LLC



Copyright ©2024 Genuin Biotechnologies LLC



Western blotting analysis using anti-Exportin 5 antibody (Cat#2741). Total cell lysates (30 μ g) from various cell lines were loaded and separated by SDS-PAGE. The blot was incubated with anti-Exportin 5 antibody (Cat#2741, 1:5,000) and HRP-conjugated goat anti-rabbit secondary antibody (Cat#201, 1:20,000) respectively. Image was developed using FeQ™ ECL Substrate Kit (Cat#226).