

Catalog #: 2861

Aliases

DENR; Density Regulated Re-Initiation And Release Factor; SMAP-3; DRP1; DRP; Density-Regulated Protein; Smooth Muscle Cell Associated Protein-3; Smooth Muscle Cell-Associated Protein 3; Protein DRP1

Background

Gene Name: DENR

NCBI Gene Entry: [8562](#)

UniProt Entry: [O43583](#)

Application Information

Molecular Weight: Predicted, 22 kDa; observed, 28 kDa

Clonality: Mouse monoclonal antibody

Clone ID: 24GB2885

Species Reactivity: Human

Applications Tested: Western blotting (WB)

Immunogen

Recombinant protein of human DENR

Isotype

Mouse IgG2b

Storage Buffer

Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:1,000-1:5,000

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Western blotting analysis using anti-DENR antibody (Cat#2861). Total cell lysates (30 μ g) from various cell lines were loaded and separated by SDS-PAGE. The blot was incubated with anti-DENR antibody (Cat#2861, 1:5,000) and HRP-conjugated goat anti-mouse secondary antibody (Cat#101, 1:20,000) respectively. Image was developed using FeQ™ ECL Substrate Kit (Cat#226).

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM