

Catalog #: 2866

Aliases

CBFA2T3; CBFA2/RUNX1 Partner Transcriptional Co-Repressor 3; ZMYND4; MTGR2; MTG16; RUNX1T3; ETO2; Core-Binding Factor, Runt Domain, Alpha Subunit 2; Translocated To, 3; Myeloid Translocation Gene On Chromosome 16 Protein; Zinc Finger MYND Domain-Containing Protein 4; Myeloid Translocation Gene 8 And 16b; CBFA2/RUNX1 Translocation Partner 3; MTG8-Related Protein; Protein CBFA2T3; MTG8-Related Gene; HMTG16

Background

Gene Name: CBFA2T3

NCBI Gene Entry: [863](#)

UniProt Entry: [O75081](#)

Application Information

Molecular Weight: Predicted, 71 kDa; observed, 71 kDa

Clonality: Mouse monoclonal antibody

Clone ID: 24GB2910

Species Reactivity: Human, mouse, rat

Applications Tested: Western blotting (WB)

Immunogen

Recombinant protein of human CBFA2T3

Isotype

Mouse IgG1 kappa

Storage Buffer

Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.

Storage


Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:1,000-1:5,000

Note: This product is for research use only.

Validation Data



Western blotting analysis using anti-CBFA2T3 antibody (Cat#2866). Total cell lysates (30 μ g) from various cell lines were loaded and separated by SDS-PAGE. The blot was incubated with anti-CBFA2T3 antibody (Cat#2866, 1:5,000) and HRP-conjugated goat anti-mouse secondary antibody (Cat#101, 1:20,000) respectively. Image was developed using FeQTM ECL Substrate Kit (Cat#226).

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM