

Catalog #: 2879

Aliases

LMNA; Lamin A/C; HGPS; MADA; LMN1; Mandibuloacral Dysplasia Type A; Lamin A/C-Like 1; Prelamin-A/C; Progerin; LGMD1B; CMD1A; LMNL1; PRO1; Limb Girdle Muscular Dystrophy 1B (Autosomal Dominant); Cardiomyopathy, Dilated 1A (Autosomal Dominant); Epididymis Secretory Sperm Binding Protein; Progeria 1 (Hutchinson-Gilford Type); Renal Carcinoma Antigen NY-REN-32; 70 KDa Lamin; Lamin C; CMT2B1; CDCD1; FPLD2; Lamin; CDDC; EMD2; FPLD; LDP1; LMNC; FPL; IDC; LFP

Background

Gene Name: LMNA

NCBI Gene Entry: [4000](#)

UniProt Entry: [P02545](#)

Application Information

Molecular Weight: Predicted, 74 kDa; observed, 63, 74 kDa

Clonality: Mouse monoclonal antibody

Clone ID: 24GB2970

Species Reactivity: Human, mouse, rat

Applications Tested: Western blotting (WB), flow cytometry (FCM)

Immunogen

Recombinant protein of human Lamin A/C

Isotype

Mouse IgG1 kappa

Storage Buffer

Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:2,500-1:5,000

Flow Cytometry (FCM): 1:2,000

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041


ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

Note: This product is for research use only.

Validation Data



Flow cytometric analysis of Lamin A/C expression in HeLa cells using anti-Lamin A/C antibody (Cat#2879, 1:2,000). Green, isotype control; red, Lamin A/C.


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Western blotting analysis using anti-Lamin A/C antibody (Cat#2879). Total cell lysates (30 μ g) from various cell lines were loaded and separated by SDS-PAGE. The blot was incubated with anti-Lamin A/C antibody (Cat#2879, 1:5,000) and HRP-conjugated goat anti-mouse secondary antibody (Cat#101, 1:20,000) respectively. Image was developed using FeQ™ ECL Substrate Kit (Cat#226).

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM