

Anti-Fibroblast Growth Factor Receptor 3 Recombinant Rabbit Monoclonal Antibody



Catalog #: 2964

Aliases

FGFR3; Fibroblast Growth Factor Receptor 3; JTK4; CD333; CEK2; EC 2.7.10.1; FGFR-3; ACH; Achondroplasia, Thanatophoric Dwarfism; Fibroblast Growth Factor Receptor 3-S; Hydroxyaryl-Protein Kinase; Tyrosine Kinase JTK4; CD333 Antigen; HSGFR3EX; EC 2.7.10

Background

Gene Name: FGFR3

NCBI Gene Entry: [2261](#)

UniProt Entry: [P22607](#)

Application Information

Molecular Weight: Predicted, 88 kDa; observed, 98 kDa

Clonality: Rabbit monoclonal antibody

Clone ID: 24GB3385

Species Reactivity: Human, mouse, rat

Applications Tested: Western blotting (WB), flow cytometry (FCM), immunocytochemistry (IC)

Immunogen

A synthesized peptide derived from human FGFR3

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:1,000-1:5,000

Flow Cytometry (FCM): 1:2,000

Immunocytochemistry (IC): 1:100-1:1,000

Note: This product is for research use only.

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS


SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

Anti-Fibroblast Growth Factor Receptor 3 Recombinant Rabbit Monoclonal Antibody

PAGE 2


Validation Data



Copyright ©2024 Genuin Biotechnologies LLC

Fibroblast growth factor receptor 3-Alexa Fluor[®] 647

Flow cytometric analysis of Fibroblast growth factor receptor 3 expression in C2C12 cells using anti-Fibroblast growth factor receptor 3 antibody (Cat#2964, 1:2,000). Green, isotype control; red, Fibroblast growth factor receptor 3.



Copyright ©2024 Genuin Biotechnologies LLC

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS


SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

Anti-Fibroblast Growth Factor Receptor 3 Recombinant Rabbit Monoclonal Antibody

PAGE 3

Immunocytochemical staining of C2C12 cells with anti-Fibroblast growth factor receptor 3 antibody (Cat#2964, 1:1,000) . Nuclei were stained blue with DAPI; Fibroblast growth factor receptor 3 was stained magenta with Alexa Fluor® 647. Images were taken using leica stellaris 5. Protein abundance based on laser intensity and smart gain: High. Scale bar: 20 µm.



Western blotting analysis using anti-Fibroblast growth factor receptor 3 antibody (Cat#2964). Total cell lysates (30 µg) from various cell lines were loaded and separated by SDS-PAGE. The blot was incubated with anti-Fibroblast growth factor receptor 3 antibody (Cat#2964, 1:5,000) and HRP-conjugated goat anti-rabbit secondary antibody (Cat#201, 1:20,000) respectively. Image was developed using FeQ™ ECL Substrate Kit (Cat#226).

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM