

Anti-Adenosine Deaminase RNA Specific Recombinant Rabbit Monoclonal Antibody



Catalog #: 3028

Aliases

ADAR; Adenosine Deaminase RNA Specific; ADAR1; DRADA; Double-Stranded RNA-Specific Adenosine Deaminase; G1P1; IFI4; 136 KDa Double-Stranded RNA-Binding Protein; Interferon-Inducible Protein 4; Interferon-Induced Protein 4; K88DSRBP; DSRAD; IFI-4; P136; Adenosine Deaminase Acting On RNA 1-A; DsRNA Adenosine Deaminase; DsRNA Adenosine Deaminase; EC 3.5.4.37; EC 3.5.4; AGS6; DSH

Background

Gene Name: ADAR

NCBI Gene Entry: [103](#)

UniProt Entry: [P55265](#)

Application Information

Molecular Weight: Predicted, 136 kDa; observed, 150 kDa

Clonality: Rabbit monoclonal antibody

Clone ID: 24GB3705

Species Reactivity: Human

Applications Tested: Western blotting (WB), immunocytochemistry (IC), flow cytometry (FCM)

Immunogen

A synthesized peptide derived from human ADAR1

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:1,000-1:5,000

Immunocytochemistry (IC): 1:100-1:1,000

Flow Cytometry (FCM): 1:2,000

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041


WWW.GENUINBIOTECH.COM

Anti-Adenosine Deaminase RNA Specific Recombinant Rabbit Monoclonal Antibody


PAGE 2

Note: This product is for research use only.

Validation Data



Flow cytometric analysis of Adenosine deaminase RNA specific expression in HepG2 cells using anti-Adenosine deaminase RNA specific antibody (Cat#3028, 1:2,000). Green, isotype control; red, Adenosine deaminase RNA specific.



Immunocytochemical staining of HepG2 cells with anti-Adenosine deaminase RNA specific antibody (Cat#3028, 1:1,000). Nuclei were stained blue with DAPI; Adenosine deaminase RNA

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS


SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

Anti-Adenosine Deaminase RNA Specific Recombinant Rabbit Monoclonal Antibody

PAGE 3

specific was stained magenta with Alexa Fluor® 647. Images were taken using Leica stellaris 5. Protein abundance based on laser Intensity and smart gain: High. Scale bar: 20 µm.



Copyright © 2024 Genuin Biotechnologies LLC

Western blotting analysis using anti-Adenosine deaminase RNA specific antibody (Cat#3028). Total cell lysates (30 µg) from various cell lines were loaded and separated by SDS-PAGE. The blot was incubated with anti-Adenosine deaminase RNA specific antibody (Cat#3028, 1:5,000) and HRP-conjugated goat anti-rabbit secondary antibody (Cat#201, 1:20,000) respectively. Image was developed using FeQ™ ECL Substrate Kit (Cat#226).

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM