

Anti-Sequestosome 1 Recombinant Rabbit Monoclonal Antibody



Catalog #: 3086

Aliases

SQSTM1; Sequestosome 1; P60; P62B; A170; OSIL; P62; Phosphotyrosine-Independent Ligand For The Lck SH2 Domain Of 62 KDa; EBI3-Associated Protein Of 60 KDa; Oxidative Stress Induced Like; Ubiquitin-Binding Protein P62; Autophagy Receptor P62; Sequestosome-1; EBIAP; PDB3; Phosphotyrosine Independent Ligand For The Lck SH2 Domain P62; EBI3-Associated Protein P60; Paget Disease Of Bone 3; FTDALS3; NADGP; DMRV; ZIP3; ORCA

Background

Gene Name: SQSTM1

NCBI Gene Entry: [8878](#)

UniProt Entry: [Q13501](#)

Application Information

Molecular Weight: Predicted, 48 kDa; observed, 62 kDa

Clonality: Rabbit monoclonal antibody

Clone ID: 24GB3995

Species Reactivity: Human, mouse, rat

Applications Tested: Western blotting (WB), flow cytometry (FCM), immunocytochemistry (IC)

Immunogen

A synthesized peptide derived from human p62/SQSTM1

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:2,000-1:10,000

Flow Cytometry (FCM): 1:2,000

Immunocytochemistry (IC): 1:100-1:1,000

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041


WWW.GENUINBIOTECH.COM

Anti-Sequestosome 1 Recombinant Rabbit Monoclonal Antibody

PAGE 2

Note: This product is for research use only.

Validation Data



Flow cytometric analysis of p62/SQSTM1 expression in HepG2 cells using anti-p62/SQSTM1 antibody (Cat#3086, 1:2,000). Green, isotype control; red, p62/SQSTM1.


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041


ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Immunocytochemical staining of HepG2 cells with anti-p62/SQSTM1 antibody (Cat#3086, 1:1,000) . Nuclei were stained blue with DAPI; p62/SQSTM1 was stained magenta with Alexa Fluor® 647. Images were taken using Leica stellaris 5. Protein abundance based on laser Intensity and smart gain: High. Scale bar: 20 μ m.



Anti-Sequestosome 1 Recombinant Rabbit Monoclonal Antibody

PAGE 4

Western blotting analysis using anti-p62/SQSTM1 antibody (Cat#3086). Total cell lysates (30 µg) from various cell lines were loaded and separated by SDS-PAGE. The blot was incubated with anti-p62/SQSTM1 antibody (Cat#3086, 1:10,000) and HRP-conjugated goat anti-rabbit secondary antibody (Cat#201, 1:20,000) respectively. Image was developed using FeQ™ ECL Substrate Kit (Cat#226).

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM