

# Biotin-Conjugated Goat Anti-Human IgA Secondary Antibody



**Catalog #: 314**

## Application Information

Clonality: Goat Polyclonal antibody

Species Reactivity: Human

Applications Tested: ELISA (E), Western Blotting (WB), Immunohistochemistry (IHC),  
Immunocytochemistry (IC)

Conjugation: Biotin

## Immunogen

Human IgA

## Storage Buffer

Supplied in PBS (pH 7.2), and 0.05% sodium azide.

## Storage

Store at -20 °C for one year.

## Recommended Dilutions

ELISA (E): 1:5,000-1:20,000

Western Blotting (WB): 1:5,000-1:20,000

Immunohistochemistry (IHC): 1:500-1:1,000

## Validation Data

---

### SUPPORT

SUPPORT@GENUINBIOTECH.COM  
TEL: +1-540-855-7041

### ORDERS

SALES@GENUINBIOTECH.COM  
FAX: +1-540-855-7041

[WWW.GENUINBIOTECH.COM](http://WWW.GENUINBIOTECH.COM)

