

Catalog #: 3399

Aliases

CDX2; Caudal Type Homeobox 2; CDX3; Caudal-Type Homeobox Protein 2; Homeobox Protein CDX-2; CDX-3; Caudal Type Homeo Box Transcription Factor 2; Caudal Type Homeobox Transcription Factor 2; Homeobox Protein MiniCDX2; CDX2/AS

Background

Gene Name: CDX2

NCBI Gene Entry: [1045](#)

UniProt Entry: [Q99626](#)

Application Information

Molecular Weight: Predicted, 34 kDa; observed, 34 kDa

Clonality: Mouse monoclonal antibody

Clone ID: 24GB6145

Species Reactivity: Human, mouse, rat

Applications Tested: Western blotting (WB)

Immunogen

Recombinant protein of human CDX2

Isotype

Mouse IgG2b

Storage Buffer

Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:1,000-1:2,500

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Copyright © 2024 Genuin Biotechnologies LLC

Western blotting analysis using anti-caudal type homeobox 2 antibody (Cat#3399). Total cell lysates (30 µg) from various cell lines were loaded and separated by SDS-PAGE. The blot was incubated with anti-caudal type homeobox 2 antibody (Cat#3399, 1:2,500) and HRP-conjugated goat anti-mouse secondary antibody (Cat#101, 1:20,000) respectively. Image was developed using NaQ™ ECL Substrate Kit (Cat#716).