

Anti-MALT1 Paracaspase Mouse Monoclonal Antibody



Catalog #: 3482

Aliases

MALT1 Paracaspase; PCASP1; MLT; Mucosa-Associated Lymphoid Tissue Lymphoma Translocation Protein 1; Mucosa Associated Lymphoid Tissue Lymphoma Translocation Gene 1; MALT1 Protease; MALT Lymphoma-Associated Translocation; Caspase-Like Protein; Paracaspase 1; Paracaspase-1; Paracaspase; EC 3.4.22.-; IMD12; MLT1

Background

Gene Name: MALT1

NCBI Gene Entry: [10892](#)

UniProt Entry: [Q9UDY8](#)

Application Information

Molecular Weight: Predicted, 92 kDa; observed, 75 kDa

Clonality: Mouse monoclonal antibody

Clone ID: 24GB5835

Species Reactivity: Human, mouse, rat

Applications Tested: Western blotting (WB)

Immunogen

Recombinant protein of human MALT1

Isotype

Mouse IgG2a

Storage Buffer

Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:2,500-1:5,000

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Western blotting analysis using anti-MALT1 antibody (Cat#3482). Total cell lysates (30 μ g) from various cell lines were loaded and separated by SDS-PAGE. The blot was incubated with anti-MALT1 antibody (Cat#3482, 1:5,000) and HRP-conjugated goat anti-mouse secondary antibody (Cat#101, 1:20,000) respectively. Image was developed using NaQ™ ECL Substrate Kit (Cat#716).