

FITC-Conjugated Goat Anti-Human IgG Fab Secondary Antibody



Catalog #: 369

Application Information

Clonality: Goat Polyclonal antibody

Species Reactivity: Human

Applications Tested: ELISA (E), Immunohistochemistry (IHC), Immunocytochemistry (IC), Flow Cytometry (FCM)

Conjugation: FITC

Immunogen

Human IgG Fab

Storage Buffer

Supplied in PBS (pH 7.2) containing 50% glycerol, 1% BSA and 0.05% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Immunocytochemistry (IC): 1:250-1:2,000

Flow Cytometry (FCM): 1:250-1:2,000

Validation Data

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

