

HRP-Conjugated Goat Anti-Human IgG Fab Secondary Antibody



Catalog #: 382

Application Information

Clonality: Goat Polyclonal antibody

Species Reactivity: Human

Applications Tested: ELISA (E), Western Blotting (WB), Immunohistochemistry (IHC),
Immunocytochemistry (IC)

Conjugation: HRP

Immunogen

Human IgG Fab

Storage Buffer

Supplied in PBS (pH 7.2) containing 50% glycerol, 1% BSA and 0.02% Proclin.

Storage

Store at -20 °C for one year.

Recommended Dilutions

ELISA (E): 1:5,000-1:20,000

Western Blotting (WB): 1:5,000-1:20,000

Immunohistochemistry (IHC): 1:500-1:1,000

Immunocytochemistry (IC): 1:500-1:1,000

Validation Data

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

**HRP-Conjugated Goat Anti-Human IgG Fab
Secondary Antibody**

PAGE 2



SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM