

# HRP-Conjugated Goat Anti-Human IgM mu chain Secondary Antibody



**Catalog #: 385**

## Application Information

Clonality: Goat Polyclonal antibody

Species Reactivity: Human

Applications Tested: ELISA (E), Western Blotting (WB), Immunohistochemistry (IHC), Immunocytochemistry (IC)

Conjugation: HRP

## Immunogen

Human IgM mu chain

## Storage Buffer

Supplied in PBS (pH 7.2) containing 50% glycerol, 1% BSA and 0.02% Proclin.

## Storage

Store at -20 °C for one year.

## Recommended Dilutions

ELISA (E): 1:5,000-1:20,000

Western Blotting (WB): 1:5,000-1:20,000

Immunohistochemistry (IHC): 1:500-1:1,000

Immunocytochemistry (IC): 1:500-1:1,000

## Validation Data

---

### SUPPORT

SUPPORT@GENUINBIOTECH.COM  
TEL: +1-540-855-7041

### ORDERS

SALES@GENUINBIOTECH.COM  
FAX: +1-540-855-7041

[WWW.GENUINBIOTECH.COM](http://WWW.GENUINBIOTECH.COM)

