

Anti-Farnesyl diphosphate synthase Recombinant Rabbit Monoclonal Antibody



Catalog #: 3987

Aliases

FDPS; Farnesyl Diphosphate Synthase; Farnesyl Pyrophosphate Synthetase, Dimethylallyltranstransferase, Geranyltranstransferase; (2E,6E)-Farnesyl Diphosphate Synthase; Farnesyl Pyrophosphate Synthase; FPP Synthase; FPS; Dimethylallyltranstransferase; Geranyltranstransferase; FPP Synthetase; EC 2.5.1.10; EC 2.5.1.1; KIAA1293; POROK9; FPPS

Background

Gene Name: FDPS

NCBI Gene Entry: [2224](#)

UniProt Entry: [P14324](#)

Application Information

Molecular Weight: Predicted, 48 kDa; observed, 36 kDa

Clonality: Rabbit monoclonal antibody

Clone ID: 24GB9025

Species Reactivity: Human

Applications Tested: Western blotting (WB), flow cytometry (FCM), immunocytochemistry (IC)

Immunogen

A synthesized peptide derived from human FDPS

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:1,000-1:5,000

Flow Cytometry (FCM): 1:2,000

Immunocytochemistry (IC): 1:100-1:1,000

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041


WWW.GENUINBIOTECH.COM

Anti-Farnesyl diphosphate synthase Recombinant Rabbit Monoclonal Antibody

PAGE 2

Note: This product is for research use only.


Validation Data



Copyright ©2025 Genuin Biotechnologies LLC

Farnesyl diphosphate synthase-Alexa Fluor® 647

Flow cytometric analysis of Farnesyl diphosphate synthase expression in HepG2 cells using anti-Farnesyl diphosphate synthase antibody (Cat#3987, 1:2,000). Green, isotype control; red, Farnesyl diphosphate synthase.



Copyright ©2025 Genuin Biotechnologies LLC

Immunocytochemical staining of HepG2 cells with anti-Farnesyl diphosphate synthase antibody

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS


SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

Anti-Farnesyl diphosphate synthase Recombinant Rabbit Monoclonal Antibody

PAGE 3

(Cat#3987, 1:1,000) . Nuclei were stained blue with DAPI; Farnesyl diphosphate synthase was stained magenta with Alexa Fluor® 647. Images were taken using Leica stellaris 5. Protein abundance based on laser Intensity and smart gain: Medium. Scale bar, 20 μ m.



Western blotting analysis using anti-farnesyl diphosphate synthase antibody (Cat#3987). Total cell lysates (30 μ g) from various cell lines were loaded and separated by SDS-PAGE. The blot was incubated with anti-farnesyl diphosphate synthase antibody (Cat#3987, 1:5,000) and HRP-conjugated goat anti-rabbit secondary antibody (Cat#201, 1:20,000) respectively. Image was developed using NaQ™ ECL Substrate Kit (Cat#716).

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM