

HRP-Conjugated Goat Anti Human Secondary Antibody



Catalog #: 401

Application Information

Clonality: Polyclonal antibody

Species Reactivity: Human

Applications Tested: Western Blotting (WB), Immunocytochemistry (IC), Immunohistochemistry (IHC), ELISA

Conjugation: HRP

Immunogen

Human IgG

Storage Buffer

Supplied in PBS (pH 7.2) containing 50% glycerol, 1% BSA and 0.02% Proclin.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:5,000-1:20,000

Immunohistochemistry (IHC): 1:500-1:1,000

Immunocytochemistry (IC): 1:500-1:1,000

Validation Data

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

