

# Texas Red-Conjugated Goat Anti-Human IgG F(ab')<sub>2</sub> Secondary Antibody



**Catalog #: 410**

## Application Information

Clonality: Goat Polyclonal antibody

Species Reactivity: Human

Applications Tested: ELISA (E), Immunohistochemistry (IHC), Immunocytochemistry (IC), Flow Cytometry (FCM)

Conjugation: Texas Red

## Immunogen

Human IgG F(ab')<sub>2</sub>

## Storage Buffer

Supplied in PBS (pH 7.2) containing 50% glycerol, 1% BSA and 0.02% sodium azide.

## Storage

Store at -20 °C for one year.

## Recommended Dilutions

Immunohistochemistry (IHC): 1:100-1:1,000

Immunocytochemistry (IC): 1:100-1:1,000

Flow Cytometry (FCM): 1:1,000-1:4,000

## Validation Data

---

### SUPPORT

SUPPORT@GENUINBIOTECH.COM  
TEL: +1-540-855-7041

### ORDERS

SALES@GENUINBIOTECH.COM  
FAX: +1-540-855-7041

[WWW.GENUINBIOTECH.COM](http://WWW.GENUINBIOTECH.COM)

**Texas Red-Conjugated Goat Anti-Human IgG  
F(ab')<sub>2</sub> Secondary Antibody**

PAGE 2



---

**SUPPORT**

SUPPORT@GENUINBIOTECH.COM  
TEL: +1-540-855-7041

**ORDERS**

SALES@GENUINBIOTECH.COM  
FAX: +1-540-855-7041

**WWW.GENUINBIOTECH.COM**