

Anti-MAPK10 Recombinant Rabbit Monoclonal Antibody



Catalog #: 4245

Aliases

Mitogen-Activated Protein Kinase 10; JNK3; P54bSAPK; P493F12; PRKM10; Stress-Activated Protein Kinase JNK3; Stress-Activated Protein Kinase 1b; C-Jun N-Terminal Kinase 3; MAP Kinase P49 3F12; MAP Kinase 10; EC 2.7.11.24; SAPK1b; JNK3A; Stress Activated Protein Kinase Beta; JNK3 Alpha Protein Kinase; EC 2.7.11; MAPK 10; SAPK1B

Background

Gene Name: MAPK10

NCBI Gene Entry: [5602](#)

UniProt Entry: [P53779](#)

Application Information

Molecular Weight: Predicted, 53 kDa; observed, 50 kDa

Clonality: Rabbit monoclonal antibody

Clone ID: 24GB10255

Species Reactivity: Human, mouse, rat

Applications Tested: Western blotting (WB), flow cytometry (FCM), immunocytochemistry (IC)

Immunogen

A synthesized peptide derived from human JNK3

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:1,000-1:5,000

Flow Cytometry (FCM): 1:2,000

Immunocytochemistry (IC): 1:100-1:1,000

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041


WWW.GENUINBIOTECH.COM

Anti-MAPK10 Recombinant Rabbit Monoclonal Antibody


PAGE 2

Note: This product is for research use only.

Validation Data



Flow cytometric analysis of JNK3 expression in HT-1080 cells using anti-JNK3 antibody (Cat#4245, 1:2,000). Green, isotype control; red, JNK3.



SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS


SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

Anti-MAPK10 Recombinant Rabbit Monoclonal Antibody

PAGE 3

Immunocytochemical staining of HT-1080 cells with anti-JNK3 antibody (Cat#4245, 1:1,000). Nuclei were stained blue with DAPI; JNK3 was stained magenta with Alexa Fluor® 647. Images were taken using Leica stellaris 5. Protein abundance based on laser Intensity and smart gain: Medium. Scale bar, 20 μ m.



Copyright ©2025 Genuin Biotechnologies LLC

Western blotting analysis using anti-JNK3 antibody (Cat#4245). Total cell lysates (30 μ g) from various cell lines were loaded and separated by SDS-PAGE. The blot was incubated with anti-JNK3 antibody (Cat#4245, 1:5,000) and HRP-conjugated goat anti-rabbit secondary antibody (Cat#201, 1:20,000) respectively. Image was developed using NaQ™ ECL Substrate Kit (Cat#716).

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM