

Anti-Phospho-PKC alpha (T638) Recombinant Rabbit Monoclonal Antibody



Catalog #: 4286

Aliases

PRKCA; Protein Kinase C Alpha; PKCA; Protein Kinase C Alpha Type; EC 2.7.11.13; PKC-Alpha; PRKACA; PKC-A; Protein Kinase C, Alpha; Aging-Associated Gene 6; EC 2.7.11; PKCalpha; PKCI+/-; PKCα; AAG6

Background

Gene Name: PRKCA

NCBI Gene Entry: [5578](#)

UniProt Entry: [P17252](#)

Application Information

Molecular Weight: Predicted, 77 kDa; observed, 77 kDa

Clonality: Rabbit monoclonal antibody

Clone ID: 24GB10460

Species Reactivity: Human, mouse, rat

Applications Tested: Western blotting (WB), flow cytometry (FCM), immunocytochemistry (IC)

Immunogen

A synthesized peptide derived from human Phospho-PKC alpha (T638)

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:1,000-1:5,000

Flow Cytometry (FCM): 1:2,000

Immunocytochemistry (IC): 1:100-1:1,000

Note: This product is for research use only.

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS


SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

Anti-Phospho-PKC alpha (T638) Recombinant Rabbit Monoclonal Antibody

PAGE 2


Validation Data



Copyright ©2024 Genuin Biotechnologies LLC

Phospho-PKC alpha (T638)-Alexa Fluor® 647

Flow cytometric analysis of phospho-PKC alpha (T638) expression in HeLa cells using anti-phospho-PKC alpha (T638) antibody (Cat#4286, 1:2,000). Green, isotype control; red, phospho-PKC alpha (T638).



Copyright ©2024 Genuin Biotechnologies LLC

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS


SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

Anti-Phospho-PKC alpha (T638) Recombinant Rabbit Monoclonal Antibody

PAGE 3

Immunocytochemical staining of Hela cells with anti-phospho-PKC alpha (T638) antibody (Cat#4286, 1:1,000). Nuclei were stained blue with DAPI; Phospho-PKC alpha (T638) was stained magenta with Alexa Fluor® 647. Images were taken using Leica stellaris 5. Protein abundance based on laser Intensity and smart gain: Medium. Scale bar, 20 µm.



Copyright ©2024 Genuin Biotechnologies LLC

Western blotting analysis using anti-phospho-PKC alpha (T638) antibody (Cat#4286). Total cell lysates (30 µg) from various cell lines were loaded and separated by SDS-PAGE. The blot was incubated with anti-phospho-PKC alpha (T638) antibody (Cat#4286, 1:5,000) and HRP-conjugated goat anti-rabbit secondary antibody (Cat#201, 1:20,000) respectively. Image was developed using NaQ™ ECL Substrate Kit (Cat#716).

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM