

Goat Anti-Human IgA Secondary Antibody



Catalog #: 438

Application Information

Clonality: Goat Polyclonal antibody

Species Reactivity: Human

Applications Tested: ELISA (E), Western Blotting (WB), Immunohistochemistry (IHC)

Immunogen

Human IgA

Storage Buffer

Supplied in PBS (pH 7.2), and 0.05% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

ELISA (E): 1:5,000-1:20,000

Western Blotting (WB): 1:5,000-1:20,000

Immunohistochemistry (IHC): 1:500-1:1,000

Validation Data

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

