

Anti-Acetyl-CoA carboxylase alpha Recombinant Rabbit Monoclonal Antibody



Catalog #: 4499

Aliases

ACACA; Acetyl-CoA Carboxylase Alpha; ACC1; ACCA; Acetyl-Coenzyme A Carboxylase Alpha; Acetyl-CoA Carboxylase 1; ACC-Alpha; HACC1; ACAC; ACACalpha; ACCalpha; Acac1; ACC; EC 6.4.1.2; ACACALPHA; ACC-ALPHA; ACCALPHA; ACACAD; ACAC1

Background

Gene Name: ACACA

NCBI Gene Entry: 31

UniProt Entry: Q13085

Application Information

Molecular Weight: Predicted, 266 kDa; observed, 266 kDa

Clonality: Rabbit monoclonal antibody

Clone ID: 24GB11525

Species Reactivity: Human, mouse, rat

Applications Tested: Western blotting (WB), flow cytometry (FCM), immunocytochemistry (IC)

Immunogen

A synthesized peptide derived from human Acetyl Coenzyme A Carboxylase

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:1,000-1:5,000

Flow Cytometry (FCM): 1:2,000

Immunocytochemistry (IC): 1:100-1:1,000

Note: This product is for research use only.

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS


SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

Anti-Acetyl-CoA carboxylase alpha Recombinant Rabbit Monoclonal Antibody

PAGE 2


Validation Data



Copyright ©2024 Genuin Biotechnologies LLC

Acetyl-CoA carboxylase alpha-Alexa Fluor® 647

Flow cytometric analysis of acetyl-CoA carboxylase alpha expression in HepG2 cells using anti-acetyl-CoA carboxylase alpha antibody (Cat#4499, 1:2,000). Green, isotype control; red, acetyl-CoA carboxylase alpha.



Copyright ©2024 Genuin Biotechnologies LLC

Immunocytochemical staining of HepG2 cells with anti-Acetyl-CoA carboxylase alpha antibody

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS


SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

Anti-Acetyl-CoA carboxylase alpha Recombinant Rabbit Monoclonal Antibody

PAGE 3

(Cat#4499, 1:1,000). Nuclei were stained blue with DAPI; Acetyl-CoA carboxylase alpha was stained magenta with Alexa Fluor® 647. Images were taken using Leica stellaris 5. Protein abundance based on laser Intensity and smart gain: Medium. Scale bar, 20 μ m.



Western blotting analysis using anti-acetyl-CoA carboxylase alpha antibody (Cat#4499). Total cell lysates (30 μ g) from various cell lines were loaded and separated by SDS-PAGE. The blot was incubated with anti-acetyl-CoA carboxylase alpha antibody (Cat#4499, 1:5,000) and HRP-conjugated goat anti-rabbit secondary antibody (Cat#201, 1:20,000) respectively. Image was developed using FeQ™ ECL Substrate Kit (Cat#226).

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM