

Anti-Legumain Recombinant Rabbit Monoclonal Antibody



Catalog #: 4506

Aliases

LGMN; Legumain; LGMN1; PRSC1; Protease, Cysteine, 1 (Legumain); Asparaginyl Endopeptidase; Protease, Cysteine 1; EC 3.4.22.34; AEP; Cysteine Protease 1

Background

Gene Name: LGMN

NCBI Gene Entry: [5641](#)

UniProt Entry: [Q99538](#)

Application Information

Molecular Weight: Predicted, 49 kDa; observed, 36,46,56 kDa

Clonality: Rabbit monoclonal antibody

Clone ID: 24GB11560

Species Reactivity: Human, mouse, rat

Applications Tested: Western blotting (WB), immunocytochemistry (IC)

Immunogen

A synthesized peptide derived from human Legumain

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:1,000-1:5,000

Immunocytochemistry (IC): 1:100-1:1,000

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041


ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Immunocytochemical staining of HeLa cells with anti-Legumain antibody (Cat#4506, 1:1,000). Nuclei were stained blue with DAPI; Legumain was stained magenta with Alexa Fluor® 647. Images were taken using Leica stellaris 5. Protein abundance based on laser Intensity and smart gain: Medium. Scale bar, 20 μ m.



Western blotting analysis using anti-legumain antibody (Cat#4506). Total cell lysates (30 μ g) from various cell lines were loaded and separated by SDS-PAGE. The blot was incubated with anti-legumain antibody (Cat#4506, 1:5,000) and HRP-conjugated goat anti-rabbit secondary antibody

Anti-Legumain Recombinant Rabbit Monoclonal Antibody

PAGE 3

(Cat#201, 1:20,000) respectively. Image was developed using FeQ™ ECL Substrate Kit (Cat#226).

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM