

Anti-Transcription Factor Dp-1 Recombinant Rabbit Monoclonal Antibody



Catalog #: 4660

Aliases

TFDP1; Transcription Factor Dp-1; DRTF1; DP1; Dp-1; DILC; Down-Regulated In Liver Cancer Stem Cells; E2F Dimerization Partner 1; DRTF1-Polypeptide 1; E2F-Related Transcription Factor

Background

Gene Name: TFDP1

NCBI Gene Entry: [7027](#)

UniProt Entry: [Q14186](#)

Application Information

Molecular Weight: Predicted, 45 kDa; observed, 49 kDa

Clonality: Rabbit monoclonal antibody

Clone ID: 24GB12310

Species Reactivity: Human, mouse, rat

Applications Tested: Western blotting (WB), flow cytometry (FCM), immunocytochemistry (IC)

Immunogen

A synthesized peptide derived from human DP1

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:1,000-1:5,000

Flow Cytometry (FCM): 1:2,000

Immunocytochemistry (IC): 1:100-1:1,000

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041


WWW.GENUINBIOTECH.COM



Copyright ©2024 Genuin Biotechnologies LLC

Transcription factor Dp-1-Alexa Fluor® 647

Flow cytometric analysis of transcription factor Dp-1 expression in C2C12 cells using anti-transcription factor Dp-1 antibody (Cat#4660, 1:2,000). Green, isotype control; red, transcription factor Dp-1.



Copyright ©2024 Genuin Biotechnologies LLC

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS


SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

Anti-Transcription Factor Dp-1 Recombinant Rabbit Monoclonal Antibody

PAGE 3

Immunocytochemical staining of C2C12 cells with anti-Transcription factor Dp-1 antibody (Cat#4660, 1:1,000). Nuclei were stained blue with DAPI; Transcription factor Dp-1 was stained magenta with Alexa Fluor® 647. Images were taken using Leica stellaris 5. Protein abundance based on laser Intensity and smart gain: Medium. Scale bar, 20 μ m.



Copyright ©2024 Genuin Biotechnologies LLC

Western blotting analysis using anti-transcription factor Dp-1 antibody (Cat#4660). Total cell lysates (30 μ g) from various cell lines were loaded and separated by SDS-PAGE. The blot was incubated with anti-transcription factor Dp-1 antibody (Cat#4660, 1:5,000) and HRP-conjugated goat anti-rabbit secondary antibody (Cat#201, 1:20,000) respectively. Image was developed using NaQ™ ECL Substrate Kit (Cat#716).

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM