

Anti-DAZ Associated Protein 1 Rabbit Monoclonal Antibody



Catalog #: 4708

Aliases

DAZAP1; DAZ Associated Protein 1; Deleted In Azoospermia-Associated Protein 1; DAZ-Associated Protein 1; MGC19907; Deleted In Azoospermia Associated Protein 1; Testicular Tissue Protein Li 50

Background

Gene Name: DAZAP1

NCBI Gene Entry: [26528](#)

UniProt Entry: [Q96EP5](#)

Application Information

Molecular Weight: Predicted, 43 kDa; observed, 43 kDa

Clonality: Rabbit monoclonal antibody

Clone ID: 24GB12550

Species Reactivity: Human, mouse, rat

Applications Tested: Western blotting (WB), flow cytometry (FCM), immunocytochemistry (IC)

Immunogen

Recombinant protein of human DAZAP1

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:1,000-1:5,000

Flow Cytometry (FCM): 1:2,000

Immunocytochemistry (IC): 1:100-1:1,000

Note: This product is for research use only.

SUPPORT


SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM


Validation Data



Copyright ©2025 Genuin Biotechnologies LLC

DAZ associated protein 1-Alexa Fluor® 647

Flow cytometric analysis of DAZ associated protein 1 expression in HepG2 cells using anti-DAZ associated protein 1 antibody (Cat#4708, 1:2,000). Green, isotype control; red, DAZ associated protein 1.



Copyright ©2025 Genuin Biotechnologies LLC

Immunocytochemical staining of HepG2 cells with anti-DAZ associated protein 1 antibody

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS


SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

Anti-DAZ Associated Protein 1 Rabbit Monoclonal Antibody

PAGE 3

(Cat#4708, 1:1,000) . Nuclei were stained blue with DAPI; DAZ associated protein 1 was stained magenta with Alexa Fluor® 647. Images were taken using Leica stellaris 5. Protein abundance based on laser Intensity and smart gain: Low. Scale bar, 20 μ m.



Western blotting analysis using anti-DAZ associated protein 1 antibody (Cat#4708). Total cell lysates (30 μ g) from various cell lines were loaded and separated by SDS-PAGE. The blot was incubated with anti-DAZ associated protein 1 antibody (Cat#4708, 1:5,000) and HRP-conjugated goat anti-rabbit secondary antibody (Cat#201, 1:20,000) respectively. Image was developed using FeQ™ ECL Substrate Kit (Cat#226).

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM