

Catalog #: 4805

Aliases

FAHD2A; Fumarylacetoacetate Hydrolase Domain Containing 2A; CGI-105; Fumarylacetoacetate Hydrolase Domain-Containing Protein 2A; Oxaloacetate Tautomerase FAHD2A, Mitochondrial; Fumarylacetoacetate Hydrolase Domain Containing 1; EC 5.3.2.2

Background

Gene Name: FAHD2A

NCBI Gene Entry: [51011](#)

UniProt Entry: [Q96GK7](#)

Application Information

Molecular Weight: Predicted, 35 kDa; observed, 31 kDa

Clonality: Mouse monoclonal antibody

Clone ID: 24GB13030

Species Reactivity: Human

Applications Tested: Western blotting (WB)

Immunogen

Recombinant protein of human FAHD2A

Isotype

Mouse IgG2b

Storage Buffer

Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:1,000-1:2,500

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Western blotting analysis using anti-FAHD2A antibody (Cat#4805). Total cell lysates (30 μ g) from various cell lines were loaded and separated by SDS-PAGE. The blot was incubated with anti-FAHD2A antibody (Cat#4805, 1:2,500) and HRP-conjugated goat anti-mouse secondary antibody (Cat#101, 1:20,000) respectively. Image was developed using FeQTM ECL Substrate Kit (Cat#226).

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM