

HRP-Conjugated Goat Anti Rabbit IgG Fc specific Secondary Antibody



Catalog #: 501

Application Information

Clonality: Polyclonal antibody

Species Reactivity: Rabbit

Applications Tested: Western Blotting (WB), Immunohistochemistry (IHC)

Conjugation: HRP

Immunogen

Rabbit IgG Fc fragment

Isotype

Goat IgG

Storage Buffer

Supplied in PBS (pH 7.4) containing 50% glycerol, 0.05% BSA and 0.02% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:10,000-1:50,000

Immunohistochemistry (IHC): 1:500-1:1,000

Validation Data

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

HRP-Conjugated Goat Anti Rabbit IgG Fc specific Secondary Antibody

PAGE 2



SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM