

Catalog #: 50102

Aliases

BM28; CCNL1; CDCL1; KIAA0030; DNA replication licensing factor MCM2; Minichromosome maintenance protein 2 homolog; Nuclear protein BM28

Background

Gene Name: MCM2

NCBI Gene Entry: [4171](#)

UniProt Entry: [P49736](#)

Application Information

Molecular Weight: Predicted, 101 kDa; observed, 130 kDa

Clonality: Rabbit monoclonal antibody

Clone ID: 24GB8400

Species Reactivity: Human

Applications Tested: Western blotting (WB)

Immunogen

Recombinant protein of human MCM2

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.3) containing 50% glycerol, and 0.05% Proclin300.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:1,000

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Western blotting analysis of MCM2 expression in HeLa (A), Jurkat (B), H1792 (C) whole cell lysates.