

Anti-Smooth Muscle Actin Rabbit Monoclonal Antibody



Catalog #: 50127

Aliases

ACTSA; ACTVS; Actin, aortic smooth muscle; Alpha-actin-2; Cell growth-inhibiting gene 46 protein

Background

Gene Name: ACTA2

NCBI Gene Entry: [59](#)

UniProt Entry: [P62736](#)

Application Information

Molecular Weight: Predicted, 42 kDa; observed, 42 kDa

Clonality: Rabbit monoclonal antibody

Clone ID: 24GB8525

Species Reactivity: Human, mouse, rat

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human Smooth Muscle Actin

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS containing 50% glycerol, 0.2% BSA and 0.02% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:1,000

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Western blotting analysis of Smooth Muscle Actin expression in HEK293T (A), HepG2 (B), H1792 (C), mouse kidney (D), mouse muscle (E), rat kidney (F), rat muscle (G) whole cell lysates.