

Catalog #: 50129

Aliases

KIAA0016; Mitochondrial import receptor subunit TOM20 homolog; Mitochondrial 20 kDa outer membrane protein; Outer mitochondrial membrane receptor Tom20

Background

Gene Name: TOMM20

NCBI Gene Entry: [9804](#)

UniProt Entry: [Q15388](#)

Application Information

Molecular Weight: Predicted, 16 kDa; observed, 16 kDa

Clonality: Rabbit monoclonal antibody

Clone ID: 24GB8535

Species Reactivity: Human, mouse, rat

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human TOM20

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS containing 50% glycerol, 0.2% BSA and 0.02% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:1,000

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Western blotting analysis of TOM20 expression in K562 (A), A549 (B), THP1 (C), mouse liver (D) whole cell lysates.