

Catalog #: 50154

Aliases

ALDH8; Aldehyde dehydrogenase family 3 member B2; Aldehyde dehydrogenase 8

Background

Gene Name: ALDH3B2

NCBI Gene Entry: [222](#)

UniProt Entry: [P48448](#)

Application Information

Molecular Weight: Predicted, 42 kDa; observed, 48 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, mouse, rat, monkey

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human ALDH3B2

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.3) containing 30% glycerol, and 0.01% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:1,000

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Western blotting analysis of ALDH3B2 expression in mouse brain (A), SGC7901 (B), A549 (C), U87MG (D) whole cell lysates.