

Anti-Arginase 1 Rabbit Polyclonal Antibody



Catalog #: 50171

Aliases

Arginase-1; Liver-type arginase; Type I arginase

Background

Gene Name: ARG1

NCBI Gene Entry: [383](#)

UniProt Entry: [P05089](#)

Application Information

Molecular Weight: Predicted, 25,34,35 kDa; observed, 40 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, mouse, rat

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human Arginase 1

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.3) containing 30% glycerol, and 0.01% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:1,000

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Western blotting analysis of Arginase 1 expression in HepG2 (A), SHSY5Y (B), A549 (C) whole cell lysates.