

Catalog #: 50177

Aliases

ZAG; ZNGP1; Zinc-alpha-2-glycoprotein; Zn-alpha-2-GP; Zn-alpha-2-glycoprotein

Background

Gene Name: AZGP1

NCBI Gene Entry: [563](#)

UniProt Entry: [P25311](#)

Application Information

Molecular Weight: Predicted, 34 kDa; observed, 47 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, monkey, monkey

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human ZAG

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.3) containing 30% glycerol, and 0.01% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:1,000

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Western blotting analysis of ZAG expression in mouse kidney (A), DLD (B) whole cell lysates.