

Catalog #: 50225

Aliases

CDKN2; MLM; Cyclin-dependent kinase inhibitor 2A, isoform 4; p14ARF; p19ARF

Background

Gene Name: CDKN2A

NCBI Gene Entry: [1029](#)

UniProt Entry: [Q8N726](#)

Application Information

Molecular Weight: Predicted, 8,11,12,13,16,17 kDa; observed, 14 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human

Applications Tested: Western blotting (WB), immunoprecipitation (IP)

Immunogen

A synthesized peptide derived from human p14 ARF

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.3) containing 30% glycerol, and 0.01% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:1000

Immunoprecipitation (IP): 1:10-1:50

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Western blotting analysis of p14 ARF expression in Saos2 (A), Molt (B), K562 (C), Jurkat (D), HL60 (E), MCF7 (F), HeLa (G) whole cell lysates.