

Catalog #: 50290

Aliases

Gasdermin-E; Inversely correlated with estrogen receptor expression 1; ICERE-1; Non-syndromic hearing impairment protein 5; DFNA5; ICERE1

Background

Gene Name: GSDME

NCBI Gene Entry: [1687](#)

UniProt Entry: [O60443](#)

Application Information

Molecular Weight: Predicted, 54 kDa; observed, 55 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, mouse, rat

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human DFNA5

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.3) containing 30% glycerol, and 0.01% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:1,000

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Western blotting analysis of DFNA5 expression in HEK293T (A) whole cell lysates.