

Catalog #: 50295

Aliases

D12S755E; DRPLA; Atrophin-1; Dentatorubral-pallidoluysian atrophy protein

Background

Gene Name: ATN1

NCBI Gene Entry: [1822](#)

UniProt Entry: [P54259](#)

Application Information

Molecular Weight: Predicted, 125 kDa; observed, 125 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, mouse, rat

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human Atrophin-1

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.3) containing 30% glycerol, and 0.01% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:1,000

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Western blotting analysis of Atrophin-1 expression in K562 (A), EOC20 (B), HeLa (C) whole cell lysates.