

Catalog #: 50308

Aliases

ENT1; Equilibrative nucleoside transporter 1; Equilibrative nitrobenzylmercaptapurine riboside-sensitive nucleoside transporter; Equilibrative NBMPR-sensitive nucleoside transporter; Nucleoside transporter, es-type; Solute carrier family 29 member 1

Background

Gene Name: SLC29A1

NCBI Gene Entry: [2030](#)

UniProt Entry: [Q99808](#)

Application Information

Molecular Weight: Predicted, 50 kDa; observed, 55,45 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, mouse, rat, bovine, dog, monkey

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human ENT1

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.3) containing 30% glycerol, and 0.01% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:1,000

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Western blotting analysis of ENT1 expression in mouse liver (A), mouse heart (B), mouse brain (C), rat heart (D), rat brain (E) whole cell lysates.