

Catalog #: 50340

Aliases

GALGT; SIAT2; Beta-1,4 N-acetylgalactosaminyltransferase 1; (N-acetylneuraminy)-galactosylglucosylceramide; GM2/GD2 synthase; GalNAc-T

Background

Gene Name: B4GALNT1

NCBI Gene Entry: [2583](#)

UniProt Entry: [Q00973](#)

Application Information

Molecular Weight: Predicted, 58 kDa; observed, 70 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, rat, monkey

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human B4GALNT1

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.3) containing 30% glycerol, and 0.01% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:1,000

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Western blotting analysis of B4GALNT1 expression in HeLa (A), C6 (B) whole cell lysates.