

Catalog #: 50358

Aliases

Probable G-protein coupled receptor 34

Background

Gene Name: GPR34

NCBI Gene Entry: [2857](#)

UniProt Entry: [Q9UPC5](#)

Application Information

Molecular Weight: Predicted, 43 kDa; observed, 44 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, mouse, rat, monkey, pig

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human GPR34

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.3) containing 30% glycerol, and 0.01% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:1,000

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Western blotting analysis of GPR34 expression in mouse brain (A), mouse kidney (B) whole cell lysates.