

Catalog #: 50387

Aliases

HOX2C; Homeobox protein Hox-B7; Homeobox protein HHO.C1; Homeobox protein Hox-2C

Background

Gene Name: HOXB7

NCBI Gene Entry: [3217](#)

UniProt Entry: [P09629](#)

Application Information

Molecular Weight: Predicted, 24 kDa; observed, 28 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, mouse, rat, bovine

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human HOXB7

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.3) containing 30% glycerol, and 0.01% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:1,000

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Western blotting analysis of HOXB7 expression in A375 (A), mouse spleen (B), rat spleen (C) whole cell lysates.