

Catalog #: 50404

Aliases

BHLHB24; ID; DNA-binding protein inhibitor ID-1; Class B basic helix-loop-helix protein 24;
bHLHb24; Inhibitor of DNA binding 1

Background

Gene Name: ID1

NCBI Gene Entry: [3397](#)

UniProt Entry: [P41134](#)

Application Information

Molecular Weight: Predicted, 16 kDa; observed, 16 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, mouse, rat, rabbit

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human ID1

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.3) containing 30% glycerol, and 0.01% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:1,000

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Western blotting analysis of ID1 expression in mouse testis (A), rat kidney (B) whole cell lysates.