

Catalog #: 50445

Aliases

CD158B1; NKAT6; Killer cell immunoglobulin-like receptor 2DL2; CD158 antigen-like family member B1; MHC class I NK cell receptor; Natural killer-associated transcript 6; NKAT-6; p58 natural killer cell receptor clone CL-43; p58 NK receptor CL-43; CD158b1

Background

Gene Name: KIR2DL2

NCBI Gene Entry: [3803](#)

UniProt Entry: [P43627](#)

Application Information

Molecular Weight: Predicted, 38 kDa; observed, 58 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, mouse, rat

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human CD158b1

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.3) containing 30% glycerol, and 0.01% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:1,000

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Western blotting analysis of CD158b1 expression in MDAMB231 (A), Jurkat (B), mouse liver (C), H9C2 (D) whole cell lysates.