

Catalog #: 50452

Aliases

KRTHB4; Keratin, type II cuticular Hb4; Keratin-84; K84; Type II hair keratin Hb4; Type-II keratin Kb24

Background

Gene Name: KRT84

NCBI Gene Entry: [3890](#)

UniProt Entry: [Q9NSB2](#)

Application Information

Molecular Weight: Predicted, 64 kDa; observed, 60 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, mouse, rat, monkey

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human Cytokeratin 84

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.3) containing 30% glycerol, and 0.01% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:1,000

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Western blotting analysis of Cytokeratin 84 expression in MCF7 (A), H1688 (B) whole cell lysates.