

Catalog #: 50489

Aliases

GST2; Microsomal glutathione S-transferase 2; Microsomal GST-2; Microsomal GST-II

Background

Gene Name: MGST2

NCBI Gene Entry: [4258](#)

UniProt Entry: [Q99735](#)

Application Information

Molecular Weight: Predicted, 16 kDa; observed, 17 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, mouse, rat, monkey

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human MGST2

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.3) containing 30% glycerol, and 0.01% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:1,000

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Western blotting analysis of MGST2 expression in MCF7 (A), SP2/0 (B), PC12 (C) whole cell lysates.