

Catalog #: 50530

Aliases

Neurogenic locus notch homolog protein 2; Notch 2; hN2

Background

Gene Name: NOTCH2

NCBI Gene Entry: [4853](#)

UniProt Entry: [Q04721](#)

Application Information

Molecular Weight: Predicted, 265 kDa; observed, 110 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, mouse, rat

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human NOTCH2

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.3) containing 30% glycerol, and 0.01% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:1,000

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Western blotting analysis of NOTCH2 expression in HEK293T (A), HepG2 (B), H1792 (C), H9C2 (D), AML12 (E) whole cell lysates.